

1986

ANNUAL QUALITY ASSURANCE PERFORMANCE REPORT

SECTION 2

BIOMATERIAL (FISH) SAMPLES

INORGANIC TRACE CONTAMINANTS SECTION

MOE
ANN
ITC
ALRK

c.1
a aa

ALRK
G. C. RONAN, DIRECTOR
Laboratory Services Branch
Ministry of the Environment

Copyright Provisions and Restrictions on Copying:

This Ontario Ministry of the Environment work is protected by Crown copyright (unless otherwise indicated), which is held by the Queen's Printer for Ontario. It may be reproduced for non-commercial purposes if credit is given and Crown copyright is acknowledged.

It may not be reproduced, in all or in part, for any commercial purpose except under a licence from the Queen's Printer for Ontario.

For information on reproducing Government of Ontario works, please contact ServiceOntario Publications at copyright@ontario.ca

1986

ANNUAL QUALITY ASSURANCE PERFORMANCE REPORT

SECTION 2

BIOMATERIAL (FISH) SAMPLES

INORGANIC TRACE CONTAMINANTS SECTION

D G STURGIS
and J C HIPFNER (editors)

Inorganic Trace Contaminants Section
Laboratory Services Branch
Ministry of the Environment

© 1989 HER MAJESTY THE QUEEN IN RIGHT OF ONTARIO

ALRK

INORGANIC TRACE CONTAMINANTS SECTION

SUMMARY

I. Introduction

The Inorganic Trace Contaminants Section of the Ministry of the Environment, Laboratory Services Branch is responsible for the analysis of a wide variety of sample types for metals and non-metals. The use of sensitive instrumentation and methodologies appropriate to the sample matrix, combined with quality assurance programs, ensures that the Section is able to maintain a high standard of analytical performance. This performance is monitored through regular internal quality control and assurance programs as well as participation in interlaboratory round-robins. This QA report summarizes the methodologies used for analysis of these samples and the supporting internal quality assurance data.

This report is assembled in sections that reflect the analyses performed on different sample matrices in support of the programs of the Ministry of the Environment. Coincidentally, these divisions also reflect the supervisory responsibilities within the Section.

II. Quality Control and Assurance

The objectives of the quality control and assurance programs are to ensure that all of the components of the analytical process are under control and to ensure immediate detection and correction of unacceptable analytical performance. The program monitors all of the reagents, instrumentation, calibration and recovery components of the analytical system.

A. Quality Control

Quality control of the analytical process takes place at the instrument level and is intended to ensure that the instrumentation is operating according to established criteria. This control function ensures that instrument calibration, standardization, slope and intercept, and instrumental drift meet these criteria.

B. Quality Assurance

Quality assurance of the analytical process takes place after the results have been generated and is intended to ensure that the analytical protocols of sample preparation and digestion have been carried out correctly. This control function ensures that reagent blanks, digested standards, sample duplicates and recovery materials meet established response criteria.

III. Report Format

The report consists of one page method summaries and one page data summaries of blanks, between-run controls and within-run duplicates in formats that are common to all of the parameter/matrix combinations. The method summaries give a brief outline of the sample preparation and measurement procedures. The data summaries consist of annual mean values with standard deviations.

For the within-run duplicates, the data set is subdivided into ranges approximating 0 to 20 %, 20 to 50 % and 50 to 100% of the analytical range. All results for duplicates reported to the data base that are "<" or that have been diluted into the range are excluded from the statistical analysis.

The standard deviations for blanks and between-run controls are calculated using formula I. Formula II is used for the calculations for within-run duplicates.

$$sd = \sqrt{[(\sum x^2 - (\sum x)^2)/n/(n-1)]} \dots\dots I$$

$$sd = \sqrt{(\sum d^2/2n)} \dots\dots II$$

where : x = the individual values; n = the number of events
d = the differences between pairs of duplicates

The data is stored in a personal computer using BMB Manager II files. All data manipulations, reports generated etc, are performed using applications written in Manager Math.

TABLE OF CONTENTS

| | | PAGE NUMBER |
|------------|------------------------------------|-------------|
| SUMMARY | INTRODUCTION | i |
| | QUALITY CONTROL AND ASSURANCE | i |
| | REPORT FORMAT | ii |
| | TABLE OF CONTENTS | iii |
| SECTION 1 | APIOs SAMPLES | 1.1 |
| | Lovol Filters | 1.3 |
| | Precipitation | 1.25 |
| | Precipitation Bag leach | 1.59 |
| SECTION 2 | BIOMATERIAL SAMPLES | 2.1 |
| SECTION 3 | AIR EMISSION SAMPLES | 3.1 |
| | Hivol Filters | 3.5 |
| | Dichotomous Filters | 3.45 |
| | Dustfall Jars | 3.95 |
| SECTION 4 | LIQUID INDUSTRIAL WASTE SAMPLES | 4.1 |
| SECTION 5 | LANDFILL LEACH SAMPLES | 5.1 |
| SECTION 6 | SOIL LEACHATE SAMPLES | 6.1 |
| SECTION 7 | DRINKING AND SURFACE WATER SAMPLES | 7.1 |
| SECTION 8 | SEDIMENT AND SOIL SAMPLES | 8.1 |
| | Sediment | 8.3 |
| | Soil | 8.73 |
| SECTION 9 | MUNICIPAL WASTE SAMPLES | 9.1 |
| | Raw Sewage | 9.3 |
| | Final Effluent | 9.35 |
| | Sludge | 9.67 |
| SECTION 10 | VEGETATION SAMPLES | 10.1 |

ITC SECTION ANNUAL QA REPORT 1986

2. Biomaterials

2.1 Fish and Biota

Fish samples are collected and frozen until they can be prepared for analysis. In most cases a filet is taken, ground thoroughly and refrozen as necessary. Other smaller biotic samples are handled as is appropriate to the sample. QA samples consist of composited fish tissue or other referenced material.

TABLE 2.1

| Parameter | Collection Device | Preparation | Analysis |
|-------------------|-------------------|-------------|--------------------|
| Metals | Plastic bags | Acid digest | AAS, ICP-AES |
| Mercury | Plastic bags | Acid digest | Cold Vapour AAS |
| Hydride Metals | Plastic bags | Acid digest | AAS |

2.2 Fish and Biota Quality Assurance

Sample duplicates are prepared by taking a second aliquot from the prepared sample.

Reagent blanks are analysed with each analytical run. There are sufficient variations in the digestion acid lots that only one lot should be used in any one analytical run.

Matrix matched between-run composite samples are prepared by collecting samples in a large container. New composites are collected as the first is depleted or as the stability period expires. These composites may be spiked as necessary to provide a measureable level of analyte.

Table 2.2 gives the sample designations for the QA materials used for Fish and Biota analyses, the nature of the sample and the parameters for which it is used as a control.

TABLE 2.2

| Sample Designation | Type | Parameter |
|--------------------|--|------------|
| con 683 | Composite in-house fish | Metals, Hg |
| con 386 | Composite in-house fish | Metals, Hg |
| fc-2 | In-house whole fish homogenized | Metals, Hg |
| fc-1 | In-house whole fish homogenized, spiked | Metals, Hg |
| TORT-1 | NRC reference lobster | Metals, Hg |
| EPA 1650 | EPA freeze dried fish | Metals, Hg |

ANALYTICAL PROCEDURE
Inorganic Trace Contaminants Section

TEST NAME: Aluminum TEST CODE: ALUT SAMPLE TYPE: Fish
UNIT: Biomaterials SUPERVISOR: R. Sadana

METHOD CODE: 550AA1

REVISION NO: 84-1

DATE: January, 1984

NATURE OF LAST REVISION:

SAMPLE HANDLING:

Quantity Required- Approximately 10 g (wet)
Container- Whirl Pac (for heavy metals only)
Preservative- None
Other-

SAMPLE PREPARATION: Partial Extn.- Total Extn.-Yes % Extracted- 100

Procedure- Weigh accurately approx. 2.5 g of a snipped sample
(See Remarks) into a calibrated 25 x 200 mm Pyrex test tube and
add 7 ml HNO₃ and 1.5 ml HClO₄.

Prepare a batch of 40 tubes in this manner, including duplicates and
spikes, and load into an aluminum heating block containing 40 holes
of appropriate size.

Heat tubes at 120°C overnight (about 15 h). Increase temperature
gradually over a 1 hour period to 180°C.

Evaporate down to approx. 0.5 ml. Add double distilled water and
and dilute to 25 ml. Mix thoroughly and submit to AAS for analysis.

INTERFERENCES: None

REPORTING RESULTS: One place after decimal if <10 ug/g; whole no if >10

INSTRUMENTATION: Perkin Elmer 5000 AAS with autosampler and PET
computer interface for data message, storage and transfer to lab
central data handling computer.

Calibration Range: 0 to 20 ug/ml

Resolution: 0.01

Sensitivity: 1.0 ug/ml

Instrument Detection Limit: 0.04 ug/ml

PERFORMANCE CHARACTERISTICS:

Routine Operating Range- 1.0000 to 200.0 mg/kg

Accuracy-

Precision of Controls-

| | A | B |
|-----------|--------|--------|
| mean | 2.747 | 9.102 |
| std. dev. | 0.6167 | 0.4111 |
| R.S.D. | 22.45 | 4.52 |

Precision of Duplicates-low range mid range high range

s.d.

mean

W

T

CONTROL LIMITS:

REMARKS: Snip samples are portions of the received fillet cut with
scissors.

SUMMARY REPORT OF QUALITY CONTROL DATA

ALUMINUM IN FISH

Range = 1.0000 to 200.0 mg/kg

IN - RUN DUPLICATES

| | | | | | |
|-------|---------|-----------------|-----------------|-----------------|--------|
| Range | <1.0000 | 1.0000 to 40.00 | 40.00 to 100.00 | 100.00 to 200.0 | >200.0 |
| no. | 0 | 0 | 0 | 0 | 0 |
| s.w. | | 0.0000 | 0.0000 | 0.0000 | |
| mean | | 0.0000 | 0.0000 | 0.0000 | |

QA CONTROL SAMPLES

| SAMPLE I.D. | NO. | MEAN | STD. DEV. | R.S.D. |
|-------------|-----|-------|-----------|--------|
| fc-1 | 0 | 0.000 | 0.0000 | 0.00 |
| fc-2 | 0 | 0.000 | 0.0000 | 0.00 |
| TORT-1 | 0 | 0.000 | 0.0000 | 0.00 |
| epa-1650 | 0 | 0.000 | 0.0000 | 0.00 |

BLANKS

| BLANK I.D. | NO. | MEAN | STD. DEV. |
|------------|-----|------|-----------|
|------------|-----|------|-----------|

ANALYTICAL PROCEDURE
Inorganic Trace Contaminants Section

TEST NAME: Cadmium TEST CODE: CDUT SAMPLE TYPE: Fish
UNIT: Biomaterials SUPERVISOR: R. Sadana

METHOD CODE: 550AA1
REVISION NO: 84-1
NATURE OF LAST REVISION:

DATE: January, 1984

SAMPLE HANDLING:

Quantity Required- Approximately 10 g (wet)
Container- Whirl Pac (for heavy metals only)
Preservative- None
Other-

SAMPLE PREPARATION: Partial Extn.- Total Extn.-Yes % Extracted- 100
Procedure- Weigh accurately approx. 2.5 g of a snipped sample
(See Remarks) into a calibrated 25 x 200 mm Pyrex test tube and
add 7 ml HNO₃ and 1.5 ml HClO₄.
Prepare a batch of 40 tubes in this manner, including duplicates and
spikes, and load into an aluminum heating block containing 40 holes
of appropriate size.
Heat tubes at 120°C overnight (about 15 h). Increase temperature
gradually over a 1 hour period to 180°C.
Evaporate down to approx. 0.5 ml. Add double distilled water and
and dilute to 25 ml. Mix thoroughly and submit to AAS for analysis.

INTERFERENCES: None

REPORTING RESULTS: One place after decimal if <10 ug/g; whole no if >10
INSTRUMENTATION: Perkin Elmer 5000 AAS with autosampler and PET
computer interface for data message, storage and transfer to lab
central data handling computer.

Calibration Range: 0.0 to 0.5 ug/ml
Resolution: 0.01
Sensitivity: 0.025 ug/ml
Instrument Detection Limit: 0.004 ug/ml

PERFORMANCE CHARACTERISTICS:

Routine Operating Range- 0.04 to 5.0 mg/kg

Accuracy-

Precision of Controls-

| | A | B |
|-----------|-----------|-----------|
| mean | .055mg/kg | .195mg/kg |
| std. dev. | .022mg/kg | .043mg/kg |
| R.S.D. | 39 % | 22 % |

Precision of Duplicates-low range mid range high range

| | | | |
|------|----|----|----|
| s.d. | ND | ND | ND |
| mean | ND | ND | ND |

W .05 mg/kg

T .20 mg/kg

CONTROL LIMITS:

REMARKS: Snip samples are portions of the received fillet cut with
scissors.
Graphite furnace can be used when lower detection limits are required.

SUMMARY REPORT OF QUALITY CONTROL DATA

CD-FLAME

IN FISH

Range = .04 to 5.0 mg/kg

IN - RUN DUPLICATES

| | | | | | | | | |
|-------|------|-----|--------|------|--------|------|--------|-------|
| Range | <.04 | .04 | to1.00 | 1.00 | to2.50 | 2.50 | to5.00 | >5.00 |
| no. | 16 | | 0 | | 0 | | 0 | 0 |
| s.w. | | 0 | | 0 | | 0 | | |
| mean | | 0 | | 0 | | 0 | | |

QA CONTROL SAMPLES

| SAMPLE I.D. | NO. | MEAN | STD. DEV. | R.S.D. |
|-------------|-----|--------|-----------|--------|
| fc-1 | 16 | 0.055 | 0.0217 | 39.45 |
| fc-2 | 11 | 0.084 | 0.0156 | 18.57 |
| TORT-1 | 7 | 24.976 | 1.4792 | 5.92 |
| epa-1650 | 6 | 0.195 | 0.0426 | 21.85 |

BLANKS

| BLANK I.D. | NO. | MEAN | STD. DEV. |
|------------|-----|------|-----------|
|------------|-----|------|-----------|

DATE

ANALYTICAL PROCEDURE
Inorganic Trace Contaminants Section

TEST NAME: Chromium TEST CODE: CRUT SAMPLE TYPE: Fish
UNIT: Biomaterials SUPERVISOR: R. Sadana

METHOD CODE: 550AA1
REVISION NO: 84-1
NATURE OF LAST REVISION:

DATE: January, 1984

SAMPLE HANDLING:

Quantity Required- Approximately 10 g (wet)
Container- Whirl Pac (for heavy metals only)
Preservative- None
Other-

SAMPLE PREPARATION: Partial Extn.- Total Extn.-Yes % Extracted- 100

Procedure- Weigh accurately approx. 2.5 g of a snipped sample
(See Remarks) into a calibrated 25 x 200 mm Pyrex test tube and
add 7 ml HNO₃ and 1.5 ml HClO₄.
Prepare a batch of 40 tubes in this manner, including duplicates and
spikes, and load into an aluminum heating block containing 40 holes
of appropriate size.
Heat tubes at 120°C overnight (about 15 h). Increase temperature
gradually over a 1 hour period to 180°C.
Evaporate down to approx. 0.5 ml. Add double distilled water and
and dilute to 25 ml. Mix thoroughly and submit to AAS for analysis.

INTERFERENCES: None

REPORTING RESULTS: One place after decimal if <10 ug/g; whole no if >10

INSTRUMENTATION: Perkin Elmer 5000 AAS with autosampler and PET
computer interface for data message, storage and transfer to lab
central data handling computer.

Calibration Range: 0 to 5 ug/ml

Resolution: .01 ug/ml

Sensitivity:

Instrument Detection Limit: .01 ug/ml

PERFORMANCE CHARACTERISTICS:

Routine Operating Range- 0.1000 to 50.0 mg/kg

Accuracy-

Precision of Controls-

| | A | B |
|-----------|--------|--------|
| mean | 2.639 | 6.107 |
| std. dev. | 0.3925 | 2.6027 |
| R.S.D. | 14.87 | 42.62 |

| Precision of Duplicates-low range | mid range | high range |
|-----------------------------------|-----------|------------|
| s.d. | | |
| mean | | |

W

T

CONTROL LIMITS:

REMARKS: Snip samples are portions of the received fillet cut with
scissors.

SUMMARY REPORT OF QUALITY CONTROL DATA

CHROMIUM IN FISH

Range = 0.1000 to 50.0 mg/kg

IN - RUN DUPLICATES

| Range | <0.1000 | 0.1000 to 10.00 | 10.00 to 25.00 | 25.00 to 50.0 | >50.0 |
|-------|---------|-----------------|----------------|---------------|-------|
| no. | 0 | 0 | 0 | 0 | 0 |
| s.w. | | 0.0000 | 0.0000 | 0.0000 | |
| mean | | 0.0000 | 0.0000 | 0.0000 | |

QA CONTROL SAMPLES

| SAMPLE I.D. | NO. | MEAN | STD. DEV. | R.S.D. |
|-------------|-----|-------|-----------|--------|
| fc-1 | 16 | 2.639 | 0.3925 | 14.87 |
| fc-2 | 0 | 0.000 | 0.0000 | 0.00 |
| TORT-1 | 4 | 3.450 | 1.1121 | 32.23 |
| epa-1650 | 0 | 0.000 | 0.0000 | 0.00 |

BLANKS

| BLANK I.D. | NO. | MEAN | STD. DEV. |
|------------|-----|------|-----------|
|------------|-----|------|-----------|

ANALYTICAL PROCEDURE
Inorganic Trace Contaminants Section

TEST NAME: Copper TEST CODE: CUUT SAMPLE TYPE: Fish
UNIT: Biomaterials SUPERVISOR: R. Sadana

METHOD CODE: 550AA1
REVISION NO: 84-1
NATURE OF LAST REVISION:

DATE: January, 1984

SAMPLE HANDLING:

Quantity Required- Approximately 10 g (wet)
Container- Whirl Pac (for heavy metals only)
Preservative- None
Other-

SAMPLE PREPARATION: Partial Extn.- Total Extn.-Yes % Extracted- 100

Procedure- Weigh accurately approx. 2.5 g of a snipped sample
(See Remarks) into a calibrated 25 x 200 mm Pyrex test tube and
add 7 ml HNO₃ and 1.5 ml HClO₄.
Prepare a batch of 40 tubes in this manner, including duplicates and
spikes, and load into an aluminum heating block containing 40 holes
of appropriate size.
Heat tubes at 120°C overnight (about 15 h). Increase temperature
gradually over a 1 hour period to 180°C.
Evaporate down to approx. 0.5 ml. Add double distilled water and
and dilute to 25 ml. Mix thoroughly and submit to AAS for analysis.

INTERFERENCES: None

REPORTING RESULTS: One place after decimal if <10 ug/g; whole no if >10
INSTRUMENTATION: Perkin Elmer 5000 AAS with autosampler and PET
computer interface for data message, storage and transfer to lab
central data handling computer.

Calibration Range: 0.0 to 5.0 ug/ml
Resolution: 0.01
Sensitivity: 0.09 ug/ml
Instrument Detection Limit: 0.02 ug/ml

PERFORMANCE CHARACTERISTICS:

Routine Operating Range- 0.1000 to 50.0 mg/kg

Accuracy-

Precision of Controls-

| | A | B |
|-----------|--------|--------|
| mean | 1.615 | 5.521 |
| std. dev. | 0.2392 | 0.3376 |
| R.S.D. | 14.81 | 6.11 |

Precision of Duplicates-low range mid range high range

s.d. 0.1118
mean 0.8060

W .01 mg/kg

T .05 mg/kg

CONTROL LIMITS:

REMARKS: Snip samples are portions of the received fillet cut with
scissors.

SUMMARY REPORT OF QUALITY CONTROL DATA

COPPER

IN FISH

Range = 0.1000to 50.0 mg/kg

IN - RUN DUPLICATES

| Range | <0.1000 | 0.1000to10.00 | 10.00 to25.00 | 25.00 to50.0 | >50.0 |
|-------|---------|---------------|---------------|--------------|-------|
| no. | 0 | 15 | 0 | 0 | 0 |
| s.w. | | 0.1118 | 0.0000 | 0.0000 | |
| mean | | 0.8060 | 0.0000 | 0.0000 | |

QA CONTROL SAMPLES

| SAMPLE I.D. | NO. | MEAN | STD. DEV. | R.S.D. |
|-------------|-----|---------|-----------|--------|
| fc-1 | 30 | 1.615 | 0.2392 | 14.81 |
| fc-2 | 0 | 0.000 | 0.0000 | 0.00 |
| TORT-1 | 6 | 375.833 | 16.4489 | 4.38 |
| epa-1650 | 0 | 0.000 | 0.0000 | 0.00 |

BLANKS

| BLANK I.D. | NO. | MEAN | STD. DEV. |
|------------|-----|------|-----------|
|------------|-----|------|-----------|

ANALYTICAL PROCEDURE
Inorganic Trace Contaminants Section

TEST NAME: Lead TEST CODE: PBUT SAMPLE TYPE: Fish
UNIT: Biomaterials SUPERVISOR: R. Sadana)

METHOD CODE: 550AA1

REVISION NO: 84-1

DATE: January, 1984

NATURE OF LAST REVISION:

SAMPLE HANDLING:

Quantity Required- Approximately 10 g (wet)
Container- Whirl Pac (for heavy metals only)
Preservative- None
Other-

SAMPLE PREPARATION: Partial Extn.- Total Extn.-Yes % Extracted- 100

Procedure- Weigh accurately approx. 2.5 g of a snipped sample
(See Remarks) into a calibrated 25 x 200 mm Pyrex test tube and
add 7 ml HNO₃ and 1.5 ml HClO₄.

Prepare a batch of 40 tubes in this manner, including duplicates and
spikes, and load into an aluminum heating block containing 40 holes
of appropriate size.

Heat tubes at 120°C overnight (about 15 h). Increase temperature
gradually over a 1 hour period to 180°C.

Evaporate down to approx. 0.5 ml. Add double distilled water and
0.5 ml of 5% potassium solution to each tube and dilute to 25 ml.

Dilute sample 5X and submit to graphite furnace for analysis.

INTERFERENCES: None

REPORTING RESULTS: One place after decimal if <10 ug/g; whole no if >10

INSTRUMENTATION: Perkin Elmer 2380 or 603 AA Spectrophotometer with a
HGA400 or HGA500 furnace and AS40 autosampler.

Calibration Range: 0 to .100 ug/ml

Resolution: 0.001

Sensitivity: 0.030 ug/ml = .200 Abs.

Instrument Detection Limit: 0.001 ug/ml

PERFORMANCE CHARACTERISTICS:

Routine Operating Range- 0.05 to 1.00 mg/kg

Accuracy-

Precision of Controls-

| | A | B |
|-----------|--------|--------|
| mean | 0.408 | 0.104 |
| std. dev. | 0.1810 | 0.0641 |
| R.S.D. | 44.36 | 61.63 |

| Precision of Duplicates- | low range | mid range | high range |
|--------------------------|-----------|-----------|------------|
| s.d. | 0.0277 | 0.0939 | |
| mean | 0.0730 | 0.3180 | |

W

T

CONTROL LIMITS:

REMARKS: Snip samples are portions of the received fillet cut with
scissors.

SUMMARY REPORT OF QUALITY CONTROL DATA

LEAD

IN FISH

Range = 0.0100to 1.0 mg/kg

IN - RUN DUPLICATES

| Range | <0.0100 | 0.0100to0.20 | 0.20 to0.50 | 0.50 to1.0 | >1.0 |
|-------|---------|--------------|-------------|------------|------|
| no. | 0 | 18 | 2 | 0 | 0 |
| s.w. | | 0.0277 | 0.0939 | 0.0000 | |
| mean | | 0.0730 | 0.3180 | 0.0000 | |

QA CONTROL SAMPLES

| SAMPLE I.D. | NO. | MEAN | STD. DEV. | R.S.D. |
|-------------|-----|-------|-----------|--------|
| fc-1 | 27 | 0.408 | 0.1810 | 44.36 |
| fc-2 | 0 | 0.000 | 0.0000 | 0.00 |
| TORT-1 | 0 | 0.000 | 0.0000 | 0.00 |
| epa-1650 | 0 | 0.000 | 0.0000 | 0.00 |

BLANKS

| BLANK I.D. | NO. | MEAN | STD. DEV. |
|------------|-----|------|-----------|
|------------|-----|------|-----------|

ANALYTICAL PROCEDURE
Inorganic Trace Contaminants Section

TEST NAME: Manganese TEST CODE: MNUT SAMPLE TYPE: Fish
UNIT: Biomaterials SUPERVISOR: R. Sadana

METHOD CODE: 550AA1
REVISION NO: 84-1
NATURE OF LAST REVISION:

DATE: January, 1984

SAMPLE HANDLING:

Quantity Required- Approximately 10 g (wet)
Container- Whirl Pac (for heavy metals only)
Preservative- None
Other-

SAMPLE PREPARATION: Partial Extn.- Total Extn.-Yes % Extracted- 100
Procedure- Weigh accurately approx. 2.5 g of a snipped sample
(See Remarks) into a calibrated 25 x 200 mm Pyrex test tube and
add 7 ml HNO₃ and 1.5 ml HClO₄.
Prepare a batch of 40 tubes in this manner, including duplicates and
spikes, and load into an aluminum heating block containing 40 holes
of appropriate size.
Heat tubes at 120°C overnight (about 15 h). Increase temperature
gradually over a 1 hour period to 180°C.
Evaporate down to approx. 0.5 ml. Add double distilled water and
and dilute to 25 ml. Mix thoroughly and submit to AAS for analysis.

INTERFERENCES: None

REPORTING RESULTS: One place after decimal if <10 ug/g; whole no if >10
INSTRUMENTATION: Perkin Elmer 5000 AAS with autosampler and PET
computer interface for data message, storage and transfer to lab
central data handling computer.

Calibration Range: 0.0 to 5 ug/ml
Resolution: 0.01
Sensitivity: 0.055 ug/ml
Instrument Detection Limit: 0.04 ug/ml

PERFORMANCE CHARACTERISTICS:

Routine Operating Range- 0.1000 to 50 mg/kg

Accuracy-

Precision of Controls-

| | A | B |
|-----------|--------|--------|
| mean | 2.739 | 1.032 |
| std. dev. | 0.2937 | 0.3195 |
| R.S.D. | 10.72 | 30.96 |

| Precision of Duplicates- | low range | mid range | high range |
|--------------------------|-----------|-----------|------------|
| s.d. | 0.1556 | ND | ND |
| mean | 0.5000 | ND | ND |

W

T

CONTROL LIMITS:

REMARKS: Snip samples are portions of the received fillet cut with
scissors.

SUMMARY REPORT OF QUALITY CONTROL DATA

MANGANESE IN FISH

Range = 0.1000 to 50.0 mg/kg

IN - RUN DUPLICATES

| Range | <0.1000 | 0.1000 to 10.00 | 10.00 to 25.00 | 25.00 to 50.0 | >50.0 |
|-------|---------|-----------------|----------------|---------------|-------|
| no. | 0 | 1 | 0 | 0 | 0 |
| s.w. | | 0.1556 | 0.0000 | 0.0000 | |
| mean | | 0.5000 | 0.0000 | 0.0000 | |

QA CONTROL SAMPLES

| SAMPLE I.D. | NO. | MEAN | STD. DEV. | R.S.D. |
|-------------|-----|--------|-----------|--------|
| fc-1 | 14 | 2.739 | 0.2937 | 10.72 |
| fc-2 | 0 | 0.000 | 0.0000 | 0.00 |
| TORT-1 | 6 | 21.133 | 0.6593 | 3.12 |
| epa-1650 | 0 | 0.000 | 0.0000 | 0.00 |

BLANKS

| BLANK I.D. | NO. | MEAN | STD. DEV. |
|------------|-----|------|-----------|
|------------|-----|------|-----------|

ANALYTICAL PROCEDURE
Inorganic Trace Contaminants Section

TEST NAME: Mercury TEST CODE: HGUT SAMPLE TYPE: Fish
UNIT: Biomaterials SUPERVISOR: R. S. Sadana

METHOD CODE: TYPE: Flameless AAS
REVISION NO: Original DATE: May, 1984
NATURE OF LAST REVISION:

SAMPLE HANDLING:

Quantity Required- Approx. 50 g
Container- Glass vials
Preservative- Freezing
Other-

SAMPLE PREPARATION: Partial Extn.- Total Extn.-Yes % Extracted-
Procedure- Weigh approx. 0.250 g of ground fish tissue into
a 50 ml Folin-Wu digestion tube. Add 5 ml of acid mixture
(4:1 - H₂SO₄:HNO₃) and place the tube in an aluminum hot block
placed on a hot plate (approx 220°C). Digest for 16 hrs (overnight).
Cool, dilute to 25 ml with distilled water.
Run in batches of 50 or more.

Treat blanks and calibration standards in exactly the same manner.
Determine mercury in the entire volume. The measurement step is
automated and is based on the evolution of atomic vapour of mercury
(wavelength - 254nm) by the addition of a reducing agent.
INTERFERENCES: Water vapour; organic solvents.

REPORTING RESULTS: Two significant figures (ug/g).

INSTRUMENTATION: Automated sampler and peristaltic pump
(Technicon or Gilson). Laboratory Data Control U.V. monitor
(Pharmacia or Milton-Roy).

Calibration Range: 0 - 16 ng/ml

Resolution: 0.4 ng/ml (one division on recorder chart paper)

Sensitivity: 1.0 ng/ml reads 0.05 Absorbance (2.5 divs on chart)

Instrument Detection Limit: 0.1 ng/ml

PERFORMANCE CHARACTERISTICS:

Routine Operating Range- 0 - 2.0 ug/g

Accuracy- 110% at 0.94 ug/g (NBS Albacore tuna)

Precision of Controls-

| | A | B |
|-----------|-----------|-----------|
| mean | .216mg/kg | .785mg/kg |
| std. dev. | .015mg/kg | .037mg/kg |
| R.S.D. | 7.0 % | 4.7 % |

| Precision of Duplicates- | low range | mid range | high range |
|--------------------------|-----------|-----------|------------|
| s.d. | .027 | .050 | .034 |
| mean | .177 | .601 | 1.23 |

W 0.01 µg/g

T 0.05 µg/g

CONTROL LIMITS:

REMARKS:

- Detection Limit - 2x std. dev. of low range within-run duplicates.
- Accuracy - Ratio of mean and cert. value in ref. mat. (%).

SUMMARY REPORT OF QUALITY CONTROL DATA

MERCURY IN FISH

Range = 0.0100to 2.0 ug/g

IN - RUN DUPLICATES

| Range | <0.0100 | 0.0100to0.40 | 0.40 to1.00 | 1.00 to2.0 | >2.0 |
|-------|---------|--------------|-------------|------------|------|
| no. | 1 | 92 | 21 | 9 | 2 |
| s.w. | | 0.0272 | 0.0496 | 0.0337 | |
| mean | | 0.1770 | 0.6010 | 1.2340 | |

QA CONTROL SAMPLES

| SAMPLE I.D. | NO. | MEAN | STD. DEV. | R.S.D. |
|-------------|-----|-------|-----------|--------|
| con 683 | 127 | 0.216 | 0.0151 | 6.99 |
| con 386 | 90 | 0.785 | 0.0368 | 4.69 |

BLANKS

| BLANK I.D. | NO. | MEAN | STD. DEV. |
|------------|-----|------|-----------|
| BLK | 0 | 0 | 0 |

DATE 86/12/17

ANALYTICAL PROCEDURE
Inorganic Trace Contaminants Section

TEST NAME: Nickel TEST CODE: NIUT SAMPLE TYPE: Fish
UNIT: Biomaterials SUPERVISOR: R. Sadana

METHOD CODE: 550AA1

REVISION NO: 1

DATE: January, 1985

NATURE OF LAST REVISION:

SAMPLE HANDLING:

Quantity Required- Approximately 10 g (wet)
Container- Whirl Pac (for heavy metals only)
Preservative- None
Other-

SAMPLE PREPARATION: Partial Extn.- Total Extn.-Yes % Extracted- 100

Procedure- Weigh accurately approx. 2.5 g of a snipped sample
(See Remarks) into a calibrated 25 x 200 mm Pyrex test tube and
add 7 ml HNO₃ and 1.5 ml HClO₄.

Prepare a batch of 40 tubes in this manner, including duplicates and
spikes, and load into an aluminum heating block containing 40 holes
of appropriate size.

Heat tubes at 120°C overnight (about 15 h). Increase temperature
gradually over a 1 hour period to 180°C.

Evaporate down to approx. 0.5 ml. Add double distilled water and
and dilute to 25 ml. Mix thoroughly and submit to AAS for analysis.

INTERFERENCES: None

REPORTING RESULTS: One place after decimal if <10 ug/g; whole no if >10

INSTRUMENTATION: Perkin Elmer 5000 AAS with autosampler and PET
computer interface for data message, storage and transfer to lab
central data handling computer.

Calibration Range: 0.0 to 5.0 ug/ml

Resolution: 0.01

Sensitivity: 0.15 ug/ml

Instrument Detection Limit: 0.04 ug/ml

PERFORMANCE CHARACTERISTICS:

Routine Operating Range- 0.1000 to 50.0 mg/kg

Accuracy-

Precision of Controls-

| | A | B |
|-----------|--------|--------|
| mean | 1.257 | 0.851 |
| std. dev. | 0.4381 | 0.1558 |
| R.S.D. | 34.85 | 18.31 |

Precision of Duplicates-low range mid range high range

s.d. 0.1939

mean 0.8040

W .02 mg/kg

T .10 mg/kg

CONTROL LIMITS:

REMARKS: Snip samples are portions of the received fillet cut with
scissors.

SUMMARY REPORT OF QUALITY CONTROL DATA

NICKEL IN FISH

Range = 0.1000to 50.0 mg/kg

IN - RUN DUPLICATES

| Range | <0.1000 | 0.1000to10.00 | 10.00 to25.00 | 25.00 to50.0 | >50.0 |
|-------|---------|---------------|---------------|--------------|-------|
| no. | 0 | 4 | 0 | 0 | 0 |
| s.w. | | 0.1939 | 0.0000 | 0.0000 | |
| mean | | 0.8040 | 0.0000 | 0.0000 | |

QA CONTROL SAMPLES

| SAMPLE I.D. | NO. | MEAN | STD. DEV. | R.S.D. |
|-------------|-----|-------|-----------|--------|
| fc-1 | 26 | 1.257 | 0.4381 | 34.85 |
| fc-2 | 0 | 0.000 | 0.0000 | 0.00 |
| TORT-1 | 6 | 3.750 | 0.9022 | 24.06 |
| epa-1650 | 0 | 0.000 | 0.0000 | 0.00 |

BLANKS

| BLANK I.D. | NO. | MEAN | STD. DEV. |
|------------|-----|------|-----------|
|------------|-----|------|-----------|

ANALYTICAL PROCEDURE
Inorganic Trace Contaminants Section

TEST NAME: Zinc TEST CODE: ZNUT SAMPLE TYPE: Fish
UNIT: Biomaterials SUPERVISOR: R. Sadana

METHOD CODE: 550AA1

REVISION NO: 84-1

DATE: January, 1984

NATURE OF LAST REVISION:

SAMPLE HANDLING:

Quantity Required- Approximately 10 g (wet)
Container- Whirl Pac (for heavy metals only)
Preservative- None
Other-

SAMPLE PREPARATION: Partial Extn.- Total Extn.-Yes % Extracted- 100

Procedure- Weigh accurately approx. 2.5 g of a snipped sample
(See Remarks) into a calibrated 25 x 200 mm Pyrex test tube and
add 7 ml HNO₃ and 1.5 ml HClO₄.

Prepare a batch of 40 tubes in this manner, including duplicates and
spikes, and load into an aluminum heating block containing 40 holes
of appropriate size.

Heat tubes at 120°C overnight (about 15 h). Increase temperature
gradually over a 1 hour period to 180°C.

Evaporate down to approx. 0.5 ml. Add double distilled water and
and dilute to 25 ml. Mix thoroughly and submit to AAS for analysis.

INTERFERENCES: None

REPORTING RESULTS: One place after decimal if <10 ug/g; whole no if >10

INSTRUMENTATION: Perkin Elmer 5000 AAS with autosampler and PET
computer interface for data message, storage and transfer to lab
central data handling computer.

Calibration Range: 0.0 to 5.0 ug/ml

Resolution: 0.01

Sensitivity: 0.18 ug/ml

Instrument Detection Limit: 0.02 ug/ml

PERFORMANCE CHARACTERISTICS:

Routine Operating Range- 0.1000 to 50.0 mg/kg

Accuracy-

Precision of Controls-

| | A | B |
|-----------|--------|--------|
| mean | 34.880 | 12.883 |
| std. dev. | 3.1565 | 5.5724 |
| R.S.D. | 9.05 | 43.25 |

Precision of Duplicates-low range mid range high range

| | | | |
|------|--------|--------|--------|
| s.d. | 0.3874 | 1.4922 | 4.0415 |
|------|--------|--------|--------|

| | | | |
|------|--------|---------|---------|
| mean | 3.9470 | 17.6000 | 30.3330 |
|------|--------|---------|---------|

W .1 mg/kg

T .5 mg/kg

CONTROL LIMITS:

REMARKS: Snip samples are portions of the received fillet cut with
scissors.

SUMMARY REPORT OF QUALITY CONTROL DATA

ZINC IN FISH

Range = 0.1000 to 50.0 mg/kg

IN - RUN DUPLICATES

| Range | <0.1000 | 0.1000 to 10.00 | 10.00 to 25.00 | 25.00 to 50.0 | >50.0 |
|-------|---------|-----------------|----------------|---------------|-------|
| no. | 0 | 10 | 3 | 3 | 0 |
| s.w. | | 0.3874 | 1.4922 | 4.0415 | |
| mean | | 3.9470 | 17.6000 | 30.3330 | |

QA CONTROL SAMPLES

| SAMPLE I.D. | NO. | MEAN | STD. DEV. | R.S.D. |
|-------------|-----|---------|-----------|--------|
| fc-1 | 30 | 34.880 | 3.1565 | 9.05 |
| fc-2 | 0 | 0.000 | 0.0000 | 0.00 |
| TORT-1 | 6 | 156.167 | 3.6009 | 2.31 |
| epa-1650 | 0 | 0.000 | 0.0000 | 0.00 |

BLANKS

| BLANK I.D. | NO. | MEAN | STD. DEV. |
|------------|-----|------|-----------|
|------------|-----|------|-----------|



(6822)

MOE/ANN/ITC/ALRK

MOE/ANN/ITC/ALRK

Ontario Ministry of the En
1986 annual quality

assurance alrk

c.1 a aa